

CENTRIFUGUE TEST TUBE 50 ml - conical with base -
30x115 mm - sterile



Product Code	37967
Unit of sale	box of 500 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	No
CH-REP	No

RDM (NSIS)	1911587
CND	W050301020102
EAN/UPC	8023279379679

Description

Centrifuges test tubes in polypropylene - graduated with cap

Test tube 50 ml - conical with base (size: 30x115 mm)
Box of 500
Sterile

All test tubes have a writing area and screw cap.