

# SYRINGES 3 PIECES WITHOUT NEEDLE - 5 ml LLC



Product Code	23807
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2325139
CND	A020102020102
EAN/UPC	8032764150945

## Description

- SYRINGES WITHOUT NEEDLE
- LUER LOCK CONE
- Capacity: 5 ml

High-quality range of medical devices ensuring an optimal safety.

The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack.

Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.