

SYRINGES 3 PIECES WITHOUT NEEDLE - 100 ml CAT



Product Code	23821
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2326511
CND	A02010203
EAN/UPC	8032764158675
PARAF	924948437

Description

- SYRINGES WITHOUT NEEDLE
- CATHETER LUER CONE
- Capacity: 100 ml

High-quality range of medical devices ensuring an optimal safety.

The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack.

Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.

High capacity polypropylene syringe.
Sterile, single use, latex free syringes for ear and other purposes.