

SYRINGES 3 PIECES WITH NEEDLE 21G - 5 ml Centric Luer Cone



Product Code	23827
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2325064
CND	A020102010201
EAN/UPC	8032764150990
PARAF	915858866

Description

- SYRINGES WITH NEEDLE
- CENTRIC LUER CONE
- Capacity: 5 ml
- 21G - 0.8 x 38 mm

High-quality range of medical devices ensuring an optimal safety.

The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack.

Needle in AISI 304 stainless steel triple bevelling.

Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.