

SYRINGES 3 PIECES WITH NEEDLE 21G - 20 ml
Eccentric LC



Product Code	23834
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2326497
CND	A020102010201
EAN/UPC	8032764150617
PARAF	905324164

Description

- SYRINGES WITH NEEDLE
- ECCENTRIC LUER CONE
- Capacity: 20 ml
- 21G - 0.8 x 38 mm

High-quality range of medical devices ensuring an optimal safety.

The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack.

Needle in AISI 304 stainless steel triple bevelling.

Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.