



Product Code	23841
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2570157
CND	A0101010202
EAN/UPC	8023279238419
UMDNS	12745
PARAF	925223442

## Description

GIMA insulin pen needles 5 mm - 31G - Universal

Insulin Pen Needles are disposable, sterile needles used with insulin pens to control high blood sugar in adults and children with diabetes mellitus.

With a fine design to minimize the discomfort of injection, ultra-thin needles provide painless and a non traumatic injection due to a silicone film.

Pyrogen free.

Multilanguage box: GB, FR, IT, DE, ES, PT, GR, PL, CZ, Arabic.

Compatible with most of insulin pens available on the market (Sanofi-Aventis, Novo Nordisk, Eli Lilly, Owen Mumford).