

# T.C. GOLD KELLY SCISSORS curved - sharp - 16 cm



|               |                |
|---------------|----------------|
| Product Code  | 38521          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | No  |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 1873018        |
| CND        | L010403        |
| EAN/UPC    | 8023279385212  |
| GTIN       | 08023279385212 |
| GMDN       | 38727          |

## Description

### GOLD LINE - TUNGSTEN CARBIDE INSTRUMENTS

A top quality theatre use Tungsten Carbide Instruments made in stainless steel AISI 410 (AISI 420 for scissors) with very hard tips made in Germany to help surgeons perform precision surgeries with best results.

Highly reliable for hardness and corrosion resistance.

Instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.