

T.C. GOLD LISTER SCISSORS - 14 cm



Product Code	38530
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2780640
CND	L01040199
EAN/UPC	8023279385304
GMDN	13481

Description

GOLD LINE - TUNGSTEN CARBIDE INSTRUMENTS

A top quality theatre use Tungsten Carbide Instruments made of stainless steel AISI 410 (AISI 420 for scissors) with very hard tips made in Germany to help surgeons perform precision surgeries with best results.

Highly reliable for hardness and corrosion resistance.

Instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.