

# BD CONNECTA THREE WAY STOPCOCKS with extension tube 50 cm - 394951



Product Code	23856
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2371104
CND	A0703
EAN/UPC	8023279238563

## Description

### BD CONNECTA™ PLUS 3 - THREE WAY STOPCOCKS

BD code: 394951

3-way stopcocks with 50 cm extens. cable

BD Connecta™ Plus 3 are 3-way stopcocks for intravenous therapy and pressure monitoring, with rotating Luer lock.

Single use, sterile.

- high flow capacity
- low priming volume (dead space): 0.23 ml
- full 360° rotation
- lipid resistant: material allows the use of lipids and chemotherapeutic agents