

T.C. GOLD GERALD FORCEPS - 17 cm - rough tips



Product Code	38560
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1829278
CND	L010403
EAN/UPC	8023279385601
GMDN	62466

Description

GOLD LINE - TUNGSTEN CARBIDE INSTRUMENTS

A top quality theatre use Tungsten Carbide Instruments made in stainless steel AISI 410 (AISI 420 for scissors) with very hard tips made in Germany to help surgeons perform precision surgeries with best results.

Highly reliable for hardness and corrosion resistance.

Instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.