

# BD EMERALD SYRINGES WITHOUT NEEDLE - 2 ml

## Centric Luer Slip



Product Code	23863
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	411629
CND	A020102010202
PARAF	987475896
GTIN	30382903077275
GMDN	47017

### Description

#### BD EMERALD SYRINGES WITHOUT NEEDLES - LUER SLIP CENTRIC TIP

- Capacity: 2 ml
- BD Code: 307727
- Minimizes Waste Disposal: Has up to 30% less material than other syringes
- Flange Size and Shape: To provide stability and comfort during use
- Textured Thumb Press: To reduce slippage during administration
- Green Stopper and Bold Scale Markings: To aid in setting an accurate dose
- Clear Barrel: To allow for visualization of fluid
- Robust Retaining Ring: To help prevent accidental plunger rod pullout