

BD EMERALD SYRINGES WITH NEEDLE 21G - 5 ml

Centric Luer Slip



Product Code	23871
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	411649
CND	A020102010201
GTIN	30382903077329
GMDN	47017

Description

BD EMERALD SYRINGE - LUER SLIP CENTRIC TIP (needle disassembled)

- Capacity: 5 ml
- Needle size: 21G - 0.8x40 mm (needle disassembled)
- BD Code: 307742
- Minimizes Waste Disposal: Has up to 30% less material than other syringes
- Flange Size and Shape: To provide stability and comfort during use
- Textured Thumb Press: To reduce slippage during administration
- Green Stopper and Bold Scale Markings: To aid in setting an accurate dose
- Clear Barrel: To allow for visualization of fluid
- Robust Retaining Ring: To help prevent accidental plunger rod pullout