

PEHA 991041 HALSTED MOSQUITO ANATOMIC FORCEPS - curved - 12.5 cm



Product Code	38733
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2494920
CND	V9099
EAN/UPC	4052199218441

Description

STERILE HAEMOSTATIC FORCEPS

- Halsted Mosquito forceps anatomic curved - 12.5 cm - 991041 - box of 25

Single-use instruments made of matt brushed, high-quality stainless steel.

Made in EU, they offer similar characteristics as multi-use instruments in terms of intended performance, haptics and design.

Sterile, individually packed and ready-to-use for a variety of surgical interventions.

They allow precise handling during procedures. Peha® color lock® safety marking allows quick and reliable differentiation from multi-use instruments.

The self-adhesive traceability label on the packaging facilitates easy, fast and less error prone data entry into the patient records.

Excellent choice for health care facilities with no processing unit.

Recyclable and disposable in existing clinical waste containers after use.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, PL, GR, HU, CZ, SI, RU, Arabic.