


# PEHA 991042 KOCHER TISSUE FORCEPS - straight - 14 cm



|               |                |
|---------------|----------------|
| Product Code  | 38735          |
| Unit of sale  | box of 25 pcs  |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2494921       |
| CND        | V9099         |
| EAN/UPC    | 4052199218489 |

## Description

### STERILE HAEMOSTATIC FORCEPS

- Kocher tissue forceps straight - 14 cm - 991042- box of 25

Single-use instruments made of matt brushed, high-quality stainless steel.

Made in EU, they offer similar characteristics as multi-use instruments in terms of intended performance, haptics and design.

Sterile, individually packed and ready-to-use for a variety of surgical interventions.

They allow precise handling during procedures. Peha® color lock® safety marking allows quick and reliable differentiation from multi-use instruments.

The self-adhesive traceability label on the packaging facilitates easy, fast and less error prone data entry into the patient records.

Excellent choice for health care facilities with no processing unit.

Recyclable and disposable in existing clinical waste containers after use.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, PL, GR, HU, CZ, SI, RU, Arabic.