

PEHA 991044 MICRO-MOSQUITO ANATOMIC FORCEPS - straight - 12.5 cm



Product Code	38738
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2494923
CND	V9099
EAN/UPC	4052199219288

Description

STERILE HAEMOSTATIC FORCEPS

- Micro-Mosquito forceps anatomic straight - 12.5 cm - 991044 - box of 25

Single-use instruments made of matt brushed, high-quality stainless steel.

Made in EU, they offer similar characteristics as multi-use instruments in terms of intended performance, haptics and design.

Sterile, individually packed and ready-to-use for a variety of surgical interventions.

They allow precise handling during procedures. Peha® color lock® safety marking allows quick and reliable differentiation from multi-use instruments.

The self-adhesive traceability label on the packaging facilitates easy, fast and less error prone data entry into the patient records.

Excellent choice for health care facilities with no processing unit.

Recyclable and disposable in existing clinical waste containers after use.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, PL, GR, HU, CZ, SI, RU, Arabic.