

# PEHA 991073 ADSON ANATOMIC FORCEPS - straight - 12 cm



Product Code	38741
Unit of sale	box of 40 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2795858
CND	V9099
EAN/UPC	4052199761473

## Description

### STERILE FORCEPS

- Adson anatomic forceps straight - 12 cm - 991073 - box of 40

Single-use instruments made of matt brushed, high-quality stainless steel.

Made in EU, they offer similar characteristics as multi-use instruments in terms of intended performance, haptics and design.

Sterile, individually packed and ready-to-use for a variety of surgical interventions.

They allow precise handling during procedures. Peha® color lock® safety marking allows quick and reliable differentiation from multi-use instruments.

The self-adhesive traceability label on the packaging facilitates easy, fast and less error prone data entry into the patient records.

Excellent choice for health care facilities with no processing unit.

Recyclable and disposable in existing clinical waste containers after use.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, PL, GR, HU, CZ, SI, RU, Arabic.