

# BD PLASTIPAK SYRINGES WITHOUT NEEDLE - 100 ml Catheter Cone



Product Code	23878
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	1869866
CND	A020102020102
GTIN	30382903006053
GMDN	47017

## Description

BD PLASTIPAK SYRINGES WHITOUT NEEDLES - CATHETER CONE

- Capacity: 100 ml
- Model: Catheter cone
- BD Code: 300605