

3.0 MHz CARDIAC PROBE (PHASED ARRAY) for code 33877



Product Code	33879
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2975208
CND	Z11040182
EAN/UPC	6945121400244
UMDNS	15976

Description

- Probe for ECOCOLORDOPPLER QBit5
- Probe type: CARDIAC (PHASED ARRAY)
 - Main Frequency: 2.5 MHz
 - Frequency Range(MHz): 1.5-5.3
 - Applications: CARDIOLOGY, TCD