

MINDRAY Z50 COLOUR ULTRASOUND with 2 probe connectors



Product Code	33881
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1947294
CND	Z11040104
EAN/UPC	6936415912475
GMDN	40761

Description

MINDRAY Z50 COLOUR PORTABLE ULTRASOUND

With 15" high resolution LCD monitor, 2 probe connectors, 500GB hard disk, 4 USB ports, multiple software package, without probes.

A colour doppler system that offers a rare combination of high performance and low cost. Superior image quality, affordability, versatility and mobility in a unique portable Doppler system.

Versatility

- scanning modes: B, Colour, PW, Power, M, FreeXros M
- image display modes: B, 2B, 4B, B/C, B/M, B/PW, Triplex PW
- applications: ABD/OB/GYN/CAR/SML/URO/VAS/ORTH/EM/NERVE
- transducer variety: Convex, Linear, Micro-convex, Endocavity
- extended scanning angle with ExFov function
- high frequency linear probe for superficial scanning

Productivity

- iTouch: quick image optimization with one-click
- iZoom: instant full screen zoom with one-click
- iStation: patient database review and analysis platform
- iStorage: direct data transfer to PC
- Clipboard: easy play, preview and sharing of cine loop
- Multiple image formats: BMP, JPG, DCM, AVI, DCM, CIN

Mobility

- compact and streamlined design weighing 8.1kg
- 15" LCD monitor with 60 degree tilt function
- dual universal transducer connectors
- two level backlit keyboard for varying requirements
- can be packed in a convenient hand-carried bag for easy transportation
- 1.5 h uninterrupted scanning with rechargeable battery

Supplied with CD manual: GB, FR, IT, ES, PT, DE, RU, BR.

Technical Specifications

Display: 15" LCD high resolution 1024x768

Software display + keyboard input: GB, FR, IT, ES, PT, RU, PL, CZ, DE

Keyboard input: TR, FI, DK, NO, SE

Power: AC 100-240 V - 50/60 Hz rechargeable lithium-ion battery - on request

Size and weight: 19x41.5x37.8 cm (7.5"x16.3"x14.9"); 8.1 kg without battery