

# STERILE IRIS SCISSORS sharp/sharp - straight - 11.5 cm



Product Code	38830
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1847810
CND	K01020102
EAN/UPC	8023279388305
PARAF	963195540
GTIN	08023279388305
GMDN	58949

## Description

- Iris scissors sharp/sharp - straight - 11.5 cm

A complete range of disposable instruments made of sturdy stainless steel alternative to plastic disposable instruments. Each instrument has a blue coloured part to avoid mix with reusable instruments and is individually packed in a sterile pouch.

Sterilized by Gamma radiation.

Instructions: GB, FR, IT, ES, DE, PT, PL, NL, SE, HU, GR, Arabic.