

STERILE HAEMOSTATIC HALSTEAD-MOSQUITO FORCEPS - straight - 12.5 cm



Product Code	38875
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes

RDM (NSIS)	1847906
CND	K01020104
EAN/UPC	8023279388756
PARAF	965565359
GTIN	08023279388756

Description

- Haemostatic Halstead-Mosquito artery forceps - straight - 12.5 cm

A complete range of disposable instruments made of sturdy stainless steel alternative to plastic disposable instruments. Each instrument has a blue coloured part to avoid mix with reusable instruments and is individually packed in a sterile pouch.

Sterilized by Gamma radiation.

Instructions: GB, FR, IT, ES, DE, PT, PL, NL, SE, HU, GR, Arabic.