

STERILE MAGILL FORCEPS - 20 cm for children



Product Code	38890
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	Yes

RDM (NSIS)	1847791
CND	K01020104
EAN/UPC	8023279388909
GMDN	33089

Description

- Magill forceps - 20 cm for children

A complete range of disposable instruments made of sturdy stainless steel alternative to plastic disposable instruments. Each instrument has a blue coloured part to avoid mix with reusable instruments and is individually packed in a sterile pouch.

Sterilized by Gamma radiation.

Instructions: GB, FR, IT, ES, DE, PT, PL, NL, SE, HU, GR, Arabic.