

STERILE MAYO-HEGAR NEEDLE HOLDER - 14 cm



Product Code	38920
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes

RDM (NSIS)	1847910
CND	K01020105
EAN/UPC	8023279389203
GTIN	08023279389203
GMDN	33961

Description

- Mayo-Hegar needle holder - 14 cm

A complete range of disposable instruments made of sturdy stainless steel alternative to plastic disposable instruments. Each instrument has a blue coloured part to avoid mix with reusable instruments and is individually packed in a sterile pouch.

Sterilized by Gamma radiation.

Instructions: GB, FR, IT, ES, DE, PT, PL, NL, SE, HU, GR, Arabic.