

STERILE MAYO-HEGAR NEEDLE HOLDER - 16 cm



Product Code	38921
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1847910
CND	K01020105
EAN/UPC	8023279389210
GMDN	33961

Description

- Mayo-Hegar needle holder - 16 cm

A complete range of disposable instruments made of sturdy stainless steel alternative to plastic disposable instruments. Each instrument has a blue coloured part to avoid mix with reusable instruments and is individually packed in a sterile pouch.

Sterilized by Gamma radiation.
Instructions: GB, FR, IT, ES, DE, PT, PL, NL, SE, HU, GR, Arabic.