

# STERILE SUTURE REMOVAL PACK



Product Code	38950
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes

RDM (NSIS)	1847926
CND	K01020199
EAN/UPC	8023279389500
PARAF	977816115
GTIN	08023279389500
GMDN	33961

## Description

- SUTURE REMOVAL PACK - box of 25

### Pack includes:

Iris scissors - straight - 11.5 cm

Adson tissue forceps - toothed - 12 cm

A complete range of disposable instruments made of sturdy stainless steel alternative to plastic disposable instruments. Each instrument has a blue coloured part to avoid mix with reusable instruments and is individually packed in a sterile pouch.

Sterilized by Gamma radiation.

Instructions: GB, FR, IT, ES, DE, PT, PL, NL, SE, HU, GR, Arabic.