

AESCULAP BUTTONED TOWEL CLAMP - curved - blunt - 14 cm - BF464R



Product Code	39255
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	Medical device
------	----------------

Class	I R
-------	-----

RDM (NSIS)	2545885
------------	---------

CND	L031307
-----	---------

EAN/UPC	4038653025367
---------	---------------

Description

Aesculap surgical instruments

Buttoned towel clamp - curved - 14 cm - BF464R