


AESCULAP HEGAR-OLSEN TC DUROGRIP NEEDLE
HOLDER - straight - 16.5 cm - BM124R



Product Code	39410
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I R

RDM (NSIS)	2422205
CND	L0205
EAN/UPC	4038653028498

Description

AESCULAP TC SURGICAL INSTRUMENTS - DUROGRIP®

- Hegar-olsen TC durogrip needle holder - straight - 16.5 cm - BM124R