


# AESCULAP HEGAR-OLSEN TC DUROGRIP NEEDLE HOLDER - straight - 16.5 cm - BM124R



Product Code	39410
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	Medical device
------	----------------

Class	I R
-------	-----

RDM (NSIS)	2422205
------------	---------

CND	L0205
-----	-------

EAN/UPC	4038653028498
---------	---------------

## Description

AESCULAP TC SURGICAL INSTRUMENTS - DUROGRIP®

- Hegar-olsen TC durogrip needle holder - straight - 16.5 cm - BM124R