


AESCULAP HEGAR-OLSEN TC DUROGRIP NEEDLE HOLDER - straight - 14.5 cm - BM128R



Product Code	39411
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I R

RDM (NSIS)	2422210
CND	L0205
EAN/UPC	4038653028511

Description

AESCULAP TC SURGICAL INSTRUMENTS - DUROGRIP®

- Hegar-olsen TC durogrip needle holder - straight - 14.5 cm - BM128R