

SONY UP-D25 MD COLOUR PRINTER



Product Code	33961
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	306367
CND	Z110706
EAN/UPC	8023279339611

Description

SONY COLOUR PRINTER UP-D25 MD - 21x9.8x39.8 cm

High quality analogue A6 Dye.
Sublimation Printer for use in medical applications.
An ideal technology for making high quality, accurate photographic prints.

USB port

It incorporates two types of colour adjustment function: RGB and HSV function