

3.0 MHz CARDIAC PROBE for code 33965-6



Product Code	33963
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	385396
CND	Z119080
EAN/UPC	8023279339635
UMDNS	15976
PARAF	903630010

Description

PHENOL FREE: completely free of phenols, it assures maximum safety to the user

Probe for Chison Q9 ecolour Dopplers (codes 33965/66).

- Scanning method: cardiac (phased array)
- Frequency (MHz)
 - Main frequency: 3.0
 - Probe band width: 2.0 - 5.4
- Main application: cardiology