

LINEAR PROBE for DP-10, DP-20, DP-2200, DP-10VET



Product Code	33977
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	170156
CND	Z1104018001
EAN/UPC	6944904054452

Description

Mindray Code: 75L38EB
Scanning Method: Linear
Main Frequency: 7.5 MHz
Probe Bandwidth: 5.0 - 10.0 MHz
Applications: Abdominal, Small Parts, MSK, Vascular, Pediatric, Orthopedics