

CHISON 5.0-10 MHz TRANSRECTAL VET PROBE for ECO1/2/3/5/6



Product Code	34000
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279340006
---------	---------------

Description

Chison Trans-rectal probe, center frequency 7.5 MHz for Portable Ultrasound Machines ECO VET series.

- Chison code: L7V-A
- Scanning method: Trans-rectal
- Main Frequency: 7.5 MHz
- Frequency Range(MHz): 5,3 - 10,0
- Applications: OVINE, EQUINE AND BOVINE

- Compliant with:
- CHISON ECO 1 VET (code 33994)
 - CHISON ECO 3EXP. VET (code 33995)
 - CHISON ECO 5 VET (code 34008)