



Product Code	40030
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8032622464245

Description

The model can be dissected on the horizontal plane to show internal details. Cornea, iris, lens and vitreous body can be removed. Muscular attachments on the sclera and part of the choroid are also shown.

Technical Specifications

- Size: 16 x 17 x h 26 cm
- Weight: approx. 1.1 kg
- Parts: 6
- Magnification: 6X