



Product Code	40072
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8032622464146

Description

This 4-part model is 3X life size.

It is sectioned along the anterior plane to show internal structures, including the cardiac valves and the comparative morphology of the right and left ventricles. The right atrium is also removable.

Technical Specifications

- Size: 24 x 25 x h 28 cm
- Weight: 2.9 kg
- Parts: 4
- Magnification: 3X