



Product Code	40077
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8032622464139

Description

Enlarged 2000 times, this model shows the various types of blood cells, including red blood cells, white blood cells (lymphocytes, monocytes, neutrophils, eosinophils and basophils) and blood platelets.

Technical Specifications

- Size: 58 x 38 x h 7 cm
- Weight: 1 kg
- Parts: 1
- Magnification: 2000X