



Product Code	40100
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8032622463620

## Description

This life size model, composed of 3 parts, shows an open dissection through the median sagittal section of the pelvis. Internal structures of the female urogenital system are portrayed in great detail. Removable parts include half of the female reproductive system with the uterus opened to reveal medial and transverse sectional views.

## Technical Specifications

- Size: 25 x 16 x h 27 cm
- Weight: 1.3 kg
- Parts: 3
- Magnification: 1X