



Product Code	40157
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279401578

Description

HUMAN SKULL COLOURED - 3 parts

A detailed reproduction of a life size human skull with all structural details. The mandible is articulated, the skullcap is removable to show internal details.

All of the joints, sutures, fissures, foramina and processes are portrayed with utmost accuracy.

Made in PVC.

Technical Specifications

- Size 19 x 13.5 x h 19 cm
- Weight: 950 g
- Magnification: 1X