



| | |
|---------------|-------------------|
| Product Code | 40162 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8032622465440 |

Description

A detailed reproduction of a life size human skull with all structural details. The skullcap is removable to show internal details. The mandible is articulated and removable. All of the joints, sutures, fissure, foramina and processes are portrayed with utmost accuracy. Significant structures are numerated. All the bones are depicted in different colours, making this model a useful tool to study the correct position and relationship of all the structures.

Technical Specifications

- Size: 19.5 x 17 x h 21.5 cm
- Weight: 950 g
- Magnification: 1X