



Product Code	40163
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8032622465235

Description

A detailed reproduction of a life size human skull with all structural details. The skullcap is removable to show internal details. The mandible is articulated and removable. All of the joints, sutures, fissure, foramina and processes are portrayed with utmost accuracy. Significant structures are numerated. The left side of the skull shows the points of origin (red) and the point of insertion (blue) of the muscles.

Technical Specifications

- Size: 18 x 17 x h 21 cm
- Weight: 875 g
- Magnification: 1X