



Product Code	40165
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8032622465365

Description

A detailed reproduction of a life size human skull with all structural details. The skullcap is removable to show internal details. The mandible is articulated and removable. All of the joints, sutures, fissure, foramina and processes are portrayed with utmost accuracy. Significant structures are numerated. It depicts all structural details of the skull. It can be disassembled into skullcap and base of skull. Articulated mandible.

Technical Specifications

- Size: 11 x 11 x h 16.5 cm
- Weight: 150 g
- Magnification: 0.5X