



Product Code	34042
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	5707480043835

Description

Adult-sized training manikin that provide an optimal lifelike representation of the human anatomy. Mechanical instrument gives you instant feedback on ventilation volume, stomach inflation, chest compression depth in mm and wrong hand position.

Make intubation possible with ET tube, Guedel, Laryngeal Mask and combi-tube.

Ambu Airway man includes torso with hooks for arms, supplied in track suit jacket, training mat and carrying bag.

Size: 80 x 35 cm
Weight: 13 kg

Technical Specifications

Mouth to Mouth ventilation: NON
 Mouth to nose ventilation: NON
 Bag-valve-mask ventilation: YES
 Use of Guedel airways: YES
 Airway management: Et tube, laryngeal tube & mask and other SGAs: YES
 Jaw-thrust manoeuvre: YES
 Gastric insufflation: YES
 Ambu hygiene system: NON
 Chest compression: YES
 Max compression depth: 7 cm
 Max insufflation: 1.2 litres
 Hand position display: YES
 Infinitely adjustable chest stiffness: YES
 Pulse: manual
 Placement of defibrillation electrodes: YES
 Defibrillation: NON
 I.V. access (need IV trainer 34041): optional
 Wireless CPR software: NON
 ECG display: NON
 Carrying bag: YES
 Weight: 13 Kg