




Product Code	34045
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	5707480134939
PARAF	950834236

Description

Training manikin for education in school of all procedures for basic life support resuscitation. It can be used for realistic training of external chest compressions and the placement of AED pads. The unique hygienic system prevents cross contamination when demonstrating and training of the principals of ventilation through mouth and nose.

Guideline conform: it allows to perform CPR training following the latest international CPR guidelines.

Compact and portable: it has a low weight and compact packaging.

Easy AED training: due to the closed shaped torso, defibrillation pads can be placed to perform a realistic AED training.

Technical Specifications

Mouth to Mouth ventilation: YES
Mouth to nose ventilation: YES
Bag-valve-mask ventilation: NON
Use of Guedel airways: NON
Airway management: Et tube, laryngeal tube & mask and other SGAs: NON
Jaw-thrust manoeuvre: NON
Gastric insufflation: NON
Ambu hygiene system: YES
Insufflation bag (Ambu hygiene system): 20
Manikin face mask (Ambu hygiene system): YES
Chest compression: YES
Max compression depth: 6 cm
Max insufflation: 0.4 litres
Hand position display: NON
Infinitely adjustable chest stiffness: NON
Pulse: NON
Placement of defibrillation electrodes: YES
Defibrillation: NON
I.V. access (need IV trainer 34041): NON
Wireless CPR software: NON
ECG display: NON
Carrying bag: NON (carton box)
Weight: 0.8 Kg