

URILYZER 100 PRO URINE ANALYZER with printer



Product Code	24049
Unit of sale	1 pc
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2428523
CND	W020101020102
EAN/UPC	8023279240498
PARAF	977810617

Description

URILYZER® 100 PRO URINE ANALYZER - with printer

Urilyzer® 100 PRO is a semi-automated and compact urine test strip reader that offers a broad range of optimized features. The intuitive user interface improves data entry and allows manual entry of colour and turbidity from a predefined list or by an additional comment. Moreover, the system provides the user with a high level of customization with flexible testing and reporting options. The intelligent database-driven data management provides multiple filter options.

To be used with urine strips 24056, 24057, 24058, 24062, 24063, 24064 only.

Manual GB, FR, IT, ES. Other languages available on www.analyticon-diagnostics.com
GB, FR, IT, ES, DE, SE, DK, NL, CZ, PL, BG, RO, RU, LT, NO.
Internal software: GB, FR, IT, ES, DE, GR, NO, SE, FI, DK, HU, CZ, PL, RO, RU.

Made in Hungary.

Technical Specifications

Measurements: 50 meas./hour in normal mode, 120 meas./hour in fast mode
Memory: 3,000 patient test results and 1,000 QC Test results
Display: 3.5" QVGA Touch-screen LCD (resolution 240x320)
Printer: Low noise thermal printer with paper roll for about 280 meas
Printout: Current number measurement, Date and time, Serial No. of the meter, patient-I-D, Asterisk on pathological results, Parameters can be printed in SI, Conventional or arbitrary units
Accuracy: Meter adjustment in 5 levels according to individual requirements
Type: Reflectance photometer wavelengths: 505, 530, 620, 660 nm
Interfaces: RS232 serial port, USB A, USB B, PS2 (external keyboard, barcode reader). Ethernet port, HL7 protocols LIS2 (ASTM+), HL7, POCT1 A2
Size - weight: 208x290x80 mm (WxDxH) / 1.2 kg
Power supply: 100-240 V, 50-60 Hz