

COMBI SCREEN CONTROL SOLUTION PN - 15 ml for 24049-51



Product Code	24055
Unit of sale	box of 2 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2428577
CND	W0101050207
EAN/UPC	4260033718198
PARAF	977810629

Description

CONTROL SOLUTION - 15 ml - box of 2

Open vial stability (at 2-8°C) of 75 days or 20 determinations (whichever occurs first)