

WHEELCHAIR STRETCHER - 2 wheels



Product Code	34060
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2299193
CND	V0899
EAN/UPC	8034028012693
PARAF	903163285

Description

Specially designed for confined areas such as narrow hallways or multiple landing stairways; anywhere access is seriously limited.

Made in light aluminium with waterproof orange sheet.

Two rearstationary wheels Ø 12 cm, two back handles and four front handles, 2 adjustable 30 cm, facilitate transport on most floor types.

2 quick lock safety belts immobilize patient during transport. Locking safety latch maintains chair position during use.

Technical Specifications

Folded size: 91 x 52 x 19.5 cm
Open size: 91 x 52 x 71 cm
Wheelchair size:
- back height: 91 cm - width: 52 cm
- seat height: 48 cm - width: 52 cm
Loading: 150 Kg
Weight: 8.5 Kg