



Product Code	34066
Unit of sale	1 kit
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
--------	----

EAN/UPC	8023279340662
GTIN	08023279340662

Description

DOCTOR'S KIT 4A

STANDARD ACCESSORIES:

Multiuse bag - blue/grey
 Aeon A200 non contact infrared thermometer (GB, FR, IT, ES)
 Oxy-4 finger oximeter - blue
 Kobe palm aneroid sphygmomanometer
 Obese cuff + 1 tube bladder
 Classic dual head stetho - Y black
 Sigma-C LED otoscope - black - pouch
 GIMA glucose monitor kit GB, FR, ES, PT
 GIMA LED penlight
 Dejerine GIMA 2000 - neurological hammer
 Fast tourniquet