



Product Code	34067
Unit of sale	1 kit
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
--------	----

EAN/UPC	8023279340679
GTIN	08023279340679

## Description

### DOCTOR'S KIT 4B

#### STANDARD ACCESSORIES:

Multiuse bag - blue/grey  
 Aeon A200 non contact infrared thermometer (GB, FR, IT, ES)  
 Oxy-4 finger oximeter - blue  
 Kobe palm aneroid sphygmomanometer  
 Obese cuff + 1 tube bladder  
 Classic dual head stetho - Y black  
 Sigma-C LED otoscope - black - pouch  
 GIMA glucose monitor kit IT, DE, GR, Arabic  
 GIMA LED penlight  
 Dejerine GIMA 2000 - neurological hammer  
 Fast tourniquet