

WHEELCHAIR STRETCHER



Product Code	34070
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2210072
CND	Y122199
EAN/UPC	8023279340709
PARAF	903163309

Description

WHEELCHAIR STRETCHER

Wheelchair stretcher with varied positions. Lightweight emergency cot quickly convert to any of four chair positions for maneuvering in narrow hallways, elevators and other confined spaces. Stretcher is made of high quality aluminium alloy. Soft foam mattress and protective block side, rust proof and easy to clean. Mattress, belts, side rails, I.V. stand and fixing device are included.

Technical Specifications

Size (LxWxH): 190x57x85 cm
Packing size: 198x60x25 cm
Load: 159 kg
Net weight, Gross weight: 32 kg, 40 kg