

VACUUM MATTRESS STRETCHER PLUS



Product Code	34075
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2598606
CND	V08050103
EAN/UPC	6974442056144

Description

VACUUM MATTRESS STRETCHER PLUS

Vacuum mattresses used to transport injured patients, to prevent compression and mitigate the risk of additional injuries. High quality and durable TPU material, filled with harmless EPS particulate material and shaped to conform to the human body. They are compatible with X-rays and CT scan. Healthcare professionals can use the air pump to generate vacuum or negative pressure, facilitating easy and rapid adjustment of softness and hardness based on severity of injury. Lightweight, foldable after deflation and easy to carry. Equipped with four handles on both sides and supplied with air pump and carrying bag.

Max. load capacity: 200 kg
Unfolded size: 204x100x5 cm
Weight: 6 kg

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic